




500 ml Centrifuge Bottle

Sl	Component	Material	Qty
1	Container	PP	1
2	Closure	PP	1

CUSTOMER COPY

DATE..... 04.05.2023

PP/PP	—	96±10 gms	Ø69.3X166mm	450 ml	00	—	—	—
Material	Shrinkage	Approx Weight	Overall Size	Brimfull volume approx	Rev.no	Revised Details	Revised by	Date
	NAME	SIGNATURE	DATE	Angle of Projection	Part Name: Centrifuge Bottle 500 ml			
DRAWN	Arunopal		03.05.2023					
CHECKED	Iranna		03.05.2023					
APPROVED	Suresh		03.05.2023					



Tarsons Products Ltd.

Remarks:

This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Products Limited.

Unspecified tolerances for open dimensions

0.00 $\pm 0.05 \text{ mm}$ 0.0 $\pm 0.1 \text{ mm}$ 0 $\pm 0.5 \text{ mm}$ Angular $\pm 0.2^\circ$

Drawing Number.

544020

A4

R00

SCALE:NTS

SHEET 1 OF 1