




1000 ml Centrifuge Bottle

Sl	Component	Material
1	Container	PP
2	Closure	PP

CUSTOMER COPY

DATE..... 05.05.2023

PP/PP	—	165±20 gms	$\varnothing 96.7 \times 173 \text{ mm}$	945 ml	00
Material	Shrinkage	Approx Weight	Overall Size	Brimfull volume approx	Rev.no
	NAME	SIGNATURE	DATE	Angle of Projection	Revised Details
DRAWN	Arunopal		03.05.2023		Revised by
CHECKED	Iranna		03.05.2023		Date
APPROVED	Suresh		03.05.2023		

Part Name:

Centrifuge Bottle 1000 ml
(Long Neck)

Tarsons Products Ltd.

Remarks:

This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Products Limited.

Unspecified tolerances for open dimensions

0.00 $\pm 0.05 \text{ mm}$ 0.0 $\pm 0.1 \text{ mm}$ 0 $\pm 0.5 \text{ mm}$ Angular $\pm 0.2^\circ$

Drawing Number.

544030

A4

R00

SCALE:NTS

SHEET 1 OF 1