




1000 ml Centrifuge Bottle

Sl	Component	Material
1	Container	PP
2	Closure	PP

CUSTOMER COPY

DATE..... 05.05.2023

PP/PP	—	155±20 gms	Ø96.7X167mm	915 ml	00	—	—	—
Material	Shrinkage	Approx Weight	Overall Size	Brimfull volume approx	Rev.no	Revised Details	Revised by	Date
	NAME	SIGNATURE	DATE	Angle of Projection	Part Name: Centrifuge Bottle 1000 ml (Short Neck)			
DRAWN	Arunopal		03.05.2023					
CHECKED	Iranna		03.05.2023					
APPROVED	Suresh		03.05.2023					



Tarsons Products Ltd.

Remarks:

This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Products Limited.

Unspecified tolerances for open dimensions

0.00 ± 0.05 mm

0.0 ± 0.1 mm

0 ± 0.5 mm

Angular ± 0.2°

Drawing Number.

544040

A4

R00

SCALE:NTS

SHEET 1 OF 1