




SL No	Description	Material	QTY
1	Bottle	HDPE	1
2	Closure	PP	1

CUSTOMER COPY

DATE: 23-08-2022

Material	Shrinkage	Approx Weight	Ø93.1x214	Overall Size	Brimfull volume approx	00	customer copy drawing	Arunopal	23.08.2022
	NAME	SIGNATURE	DATE	Angle of Projection	Part Name:	Rev.no	Revised Details	Revised by	Date
DRAWN	Arunopal		23.08.2022		1000 ml Round Narrow Mouth Bottle(Amber)				
CHECKED	Iranna								
APPROVED	Suresh								



Tarsons Products Ltd.

Remarks:

This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Product Limited.

Unspecified tolerances for open dimensions

0.00 ± 0.05 mm

0.0 ± 0.1 mm

0 ± 0.5 mm

Angular ± 0.2°

Drawing Number.

581250

A4

R00

SCALE:1:1

SHEET 1 OF 1