



If in doubt please ask

All dimensions are in mm



CUSTOMER COPY

08-04-2022

DATE.....

HDPE	%	36.7 ± 5.0 gms		wxbxh			SI no	Description	Revisions	Date
MATERIAL	Shrinkage	WEIGHT		FINISH SIZE		Brimfull volume :	DO NOT SCALE DRAWING		Revision	
	NAME	SIGNATURE	DATE	TITLE:	Angle of Projection		Description 250 ml Round Wide Mouth Bottle Amber			
DRAWN	Shubham		08-04-2022							
CHK'D	Iranna				Remarks:					
APPV'D	Suresh									
DEBURR AND BREAK SHARP EDGES		Unspecified tolerances for open dimension								
		0 - 50mm	+/-0.1 mm	This is confidential document. It must not be copied wholly or in part not be disclosed to third party or not be amended without prior written consent from M/S tarsons products pvt ltd			DWG NO. 581330		A4	
		50 - 100mm	+/-0.2mm							
		100 <	+/-0.25mm							
				Angular	+/- 0.5°				SCALE:NTS	