



IF IN DOUBT PLEASE ASK
ALL DIMENSIONS ARE IN mm



SECTION A-A



CUSTOMER COPY

DATE.....20.10.2022

PP	—	12.5±1gms	—	—	00	Customer Copy	—	—
Material	Shrinkage	Approx Weight	Overall Size	Brimfull volume approx	Rev.no	Revised Details	Revised by	Date
	NAME	SIGNATURE	DATE	Angle of Projection 	Part Name: 30 ml round wide mouth bottle			
DRAWN	Arunopal		20.10.2022					
CHECKED	Iranna							
APPROVED	Suresh							

Part Name:

30 ml round wide mouth bottle

Drawing Number.

582200

A4

R00

Remarks:

Unspecified tolerances for open dimensions	
--	--

0.00	± 0.05 mm
------	---------------

0.0	± 0.1 mm
-----	--------------

0	$\pm 0.5 \text{ mm}$
---	----------------------

Angular	$\pm 0.2^\circ$
---------	-----------------

SCALE:

SHEET 1 OF 1



Tarsons Products Ltd.

This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Product Limited.