




SECTION A-A

CUSTOMER COPY

DATE 20.10.2022

PP	—	16.3±1gms	—	—	00	Customer Copy	—	—
Material	Shrinkage	Approx Weight	Overall Size	Brimfull volume approx	Rev.no	Revised Details	Revised by	Date
	NAME	SIGNATURE	DATE	Angle of Projection	Part Name: 60 ml round wide mouth bottle			
DRAWN	Arunopal		20.10.2022					
CHECKED	Iranna							
APPROVED	Suresh							



Tarsons Products Ltd.

Remarks:

This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Product Limited.

Unspecified tolerances for open dimensions	
0.00	± 0.05 mm
0.0	± 0.1 mm
0	± 0.5 mm
Angular	± 0.2°

Drawing Number.	A4
582210	R00
SCALE:	SHEET 1 OF 1