



If in doubt please ask
All dimensions are in mm



CUSTOMER COPY

DATE 08-04-2022

PP	%	36.7 ± 5.0 gms		wxbxh			Sl no	Description	Revisions	Date
MATERIAL	Shrinkage	WEIGHT		FINISH SIZE		Brimfull volume :	DO NOT SCALE DRAWING		Revision	
	NAME	SIGNATURE	DATE	TITLE:	Angle of Projection		Description 250 ml Round Wide Mouth Bottle			
DRAWN	Shubham		08-04-2022							
CHK'D	Iranna				Remarks:					
APPV'D	Suresh									
DEBURR AND BREAK SHARP EDGES		Unspecified tolerances for open dimension		This is confidential document. It must not be copied wholly or in part not be disclosed to third party or not be amended without prior written consent from M/S tarsons products pvt ltd						
		0 - 50mm	+/-0.1 mm							
		50 - 100mm	+/-0.2mm							
		100 <	+/-0.25mm							
		Angular	+/- 0.5°			SCALE:NTS		SHEET 1 OF 1		