


IF IN DOUBT PLEASE ASK
ALL DIMENSIONS ARE IN mm



500 ml Round WM Bottle(PP)			
SL	Description	MAT	QTY
1	Bottle	PP	1
2	Closure(53-415)	PP	1

CUSTOMER COPY DATE...17.07.2023.....	PP	—	67.7±1 gms	—	—	00	Customer Copy Drawing	Arunopal	17.07.2023				
	Material	Shrinkage	Approx Weight	Overall Size	Brimfull volume approx	Rev.no	Revised Details	Revised by	Date				
		NAME	SIGNATURE	DATE	Angle of Projection	Part Name: 500 ml Round Wide Mouth PP Bottle							
	DRAWN	Arunopal		17.07.2023									
	CHECKED	Iranna		17.07.2023									
	APPROVED	Suresh		17.07.2023		Drawing Number. 582240							
Remarks:					Unspecified tolerances for open dimensions								
					0.00	± 0.05 mm	A4 R00						
					0.0	± 0.1 mm							
					0	± 0.5 mm							
Tarsons Products Ltd.					Angular	± 0.2°	SCALE:NTS		SHEET 1 OF 1				
This is a confidential document. It is strictly prohibited to be disclosed to any third party without prior permission. This must not be amended without written consent from M/S Tarsons Product Limited.													