



3800 NORTH MILL ROAD
VINELAND NJ 08360
PHONE: 1-800-843-1794
FAX: 1-800-922-4361

Certificate of Compliance

We hereby certify that our Part Number **CG-3503-10**, *Thermometer, non merc, -20 -110c* has been checked to the specifications set forth by Chemglass and is free of defects.

Material: Glass, Red Organic Liquid

Dimensions: 12" (300mm) Long; Immersion Length = 3" (76mm)

Other: Student Grade Thermometer, Range -20°C to 110°C, Division 1°C, Accuracy +/- 1.5°C, +/- 2°C above 100°C
Is not traceable to NIST, No Statement of Accuracy and Traceability available

Made in China

Date of Shipment: _____

Quantity Shipped: _____

Order #: _____

Customer P.O. _____

Signed: Hany L. Custino (QA MANAGER)

Date: 6/23/21