

Screw Cap for Sample Vial, Kimble®

김블 스크류 캡, 샘플 바이알 용

• DWK Lifescience / Kimble®

Black Phenolic PTFE Faced Rubber Lined Cap

Cat. No.	Model	Cap 규격	EA Price ₩	PK Price ₩
B08-214-111	73802-8425	8-425	800	144ea 107,500
B08-214-115	73802-13415	13-415	600	125ea 70,700
B08-214-119	73802-13425	13-425	1,300	144ea 166,500
B08-214-123	73802-15415	15-415	700	125ea 77,800
B08-214-127	73802-15425	15-425	1,400	144ea 180,600
B08-214-131	73802-18400	18-400	1,900	144ea 248,600
B08-214-135	73802-20400	20-400	2,600	144ea 338,700
B08-214-139	73802-22400	22-400	2,700	144ea 362,700
B08-214-143	73802-24400	24-400	2,700	144ea 356,400
B08-214-151	73802-33430	33-430	3,900	144ea 525,100
B08-214-154	73802-38430	38-430	5,100	144ea 687,700



Black Phenolic Cemented-In Rubber Lined Cap

Cat. No.	Model	Cap 규격	PK Price ₩	CS Price ₩
B08-214-001	73800-13415	13-415	100ea 26,900	10pk 251,400
B08-214-004	73800-15415	15-415	100ea 28,300	10pk 264,000
B08-214-007	73800-18415	18-415	100ea 35,000	10pk 326,500
V05-137-519	75204G-20400	20-400	12ea 10,400	12pk 115,800
V05-137-522	75204G-22400	22-400	12ea 11,500	12pk 128,300
V05-137-525	75204G-24400	24-400	12ea 12,300	12pk 137,200
V05-137-531	75204G-28400	28-400	12ea 11,200	12pk 125,200
V05-137-537	45066B-28	28-410	-	150ea 1,031,400
V05-137-540	75204G-33400	33-400	12ea 13,900	12pk 155,100
V05-137-543	73803-33430	33-430	-	144ea 313,000
V05-137-546	75204G-38400	38-400	12ea 16,900	12pk 189,400
V05-137-549	73803-38430	38-430	-	144ea 355,200
V05-137-555	75204G-43400	43-400	12ea 18,200	12pk 204,000
V05-137-558	75204G-45400	45-400	12ea 23,600	12pk 264,200
V05-137-561	75204G-48400	48-400	12ea 25,700	12pk 288,000
V05-137-564	75204G-53400	53-400	12ea 30,400	12pk 340,100
V05-137-567	75204G-58400	58-400	12ea 34,300	12pk 384,200
V05-137-570	75204G-63400	63-400	12ea 42,300	12pk 474,300
V05-137-573	75204G-70400	70-400	12ea 53,400	12pk 598,200
V05-137-576	75204G-89400	89-400	12ea 69,900	12pk 783,500




• Autoclavable

• Identified in ASTM Specification E892, Class A Requirements

White Urea PTFE Faced Foam Backed Rubber Lined Cap

Cat. No.	Model	Cap 규격	PK Price ₩	CS Price ₩
V05-137-601	73802U-8425	8-425	144ea 114,200	3pk 320,100
V05-137-605	73802U-13425	13-425	144ea 82,700	3pk 231,800
V05-137-609	73802U-15425	15-425	144ea 113,600	3pk 318,300
V05-137-613	73802U-18400	18-400	144ea 117,100	3pk 328,200
V05-137-617	73802U-20400	20-400	144ea 146,900	3pk 411,700
V05-137-621	73802U-24400	24-400	144ea 372,700	3pk 1,044,900



Natural PP Linerless Cap

플라스틱 캡

Cat. No.	Model	GPI	PK Price ₩	CS Price ₩
B08-214-051	73805B-13415	13-415	100ea 18,900	10pk 176,700
B08-214-054	73805B-15415	15-415	100ea 19,500	10pk 181,700
B08-214-057	73805B-18415	18-415	100ea 24,900	5pk 116,000

• 재질 : Natural PP, -10 ~ 135°C 내열, Autoclavable

• DWK Lifescience / Kimble®

