

PTFE Screw Joint for Tubing Fitting PTFE 테프론 튜빙 연결 피팅

- 재질 : PTFE, -200 ℃ to 250℃ 내열
- 사용 압력 : Max 5 bar
- 사용 가능 튜빙 및 튜브 재질 : PFA, PTFE, FEP or Glass, Steel etc.
- 테프론 재질로 최상의 내화학성을 가짐
- 시료는 단지 PTFE만 접촉 되어짐
- BOLA

Union, 일자형

Cat. No.	Model	Thread	Tubing	Bore Dia	Total Length	EA Price ₩
C20-460-201	D503-02	GL 14	od 4 mm	Ø 6 mm	49mm	41,600
C20-460-205	D503-04		od 6 mm			40,500
C20-460-209	D503-06		od 6.35(1/4") mm			44,500
C20-460-213	D503-08	GL 18	od 8 mm	Ø 8 mm	54mm	54,700
C20-460-217	D503-12		od 10 mm			56,400
C20-460-221	D503-14	GL 28	od 12 mm	Ø 14 mm	58mm	147,100
C20-460-225	D503-16		od 14 mm			144,700
C20-460-229	D503-18		od 16 mm			155,900

Union Elbow, L자형

Cat. No.	Model	Thread	Tubing	Bore Dia	Dimension	EA Price ₩
C20-460-251	D504-02	GL 14	od 4 mm	Ø 4 mm	39 x 39mm	78,700
C20-460-255	D504-04		od 6 mm			77,100
C20-460-259	D504-08	GL 18	od 8 mm	Ø 8 mm	43 x 43mm	88,600
C20-460-263	D504-12		od 10 mm			90,400
C20-460-267	D504-14	GL 28	od 12 mm	Ø 14 mm	54 x 54mm	144,900
C20-460-271	D504-16		od 14 mm			143,700
C20-460-275	D504-18		od 16 mm			161,600

Reducing Union 축소형

Cat. No.	Model	Thread ①	Tubing ①	Bore Dia	Thread ②	Tubing ②	EA Price ₩
C20-460-301	D526-02	GL 14	od 6 mm	Ø 6 mm	GL 14	od 4 mm	58,300
C20-460-305	D526-04	GL 18	od 8 mm			od 6 mm	65,100
C20-460-309	D526-10		od 6 mm			63,400	
C20-460-313	D526-06		od 4 mm			65,700	
C20-460-317	D526-12		od 6 mm			63,700	
C20-460-321	D526-14		Ø 8 mm	GL 18	od 8 mm	63,000	
C20-460-325	D526-26	GL 28	od 12 mm	Ø 6 mm	GL14	od 4 mm	152,700
C20-460-329	D526-32			Ø 10 mm	GL18	od 6 mm	140,700
C20-460-333	D526-38					od 8 mm	156,500
C20-460-337	D526-18					od 10 mm	124,500
C20-460-341	D526-28		od 14 mm	Ø 6 mm	GL14	od 4 mm	150,500
C20-460-345	D526-34			Ø 10 mm	GL18	od 6 mm	148,700
C20-460-349	D526-20		od 16 mm			od 10 mm	122,300
C20-460-353	D526-30		od 16 mm	Ø 6 mm	GL 14	od 4 mm	159,500
C20-460-357	D526-36					od 6 mm	158,000
C20-460-361	D526-22					Ø 10 mm	GL 18

T-Type

Cat. No.	Model	Thread	Tubing	Bore Dia	Dimension	EA Price ₩
C20-460-375	D505-02	GL 14	od 4 mm	Ø 4 mm	56 x 39mm	85,700
C20-460-378	D505-04		od 6 mm			84,100
C20-460-381	D505-06		od 6.35(1/4") mm			89,900
C20-460-384	D505-08	GL 18	od 8 mm	Ø 8 mm	60 x 43mm	108,700
C20-460-387	D505-12		od 10 mm			111,200
C20-460-390	D505-14	GL 28	od 12 mm	Ø 14 mm	71 x 54mm	156,400
C20-460-393	D505-16		od 14 mm			155,300
C20-460-396	D505-18		od 16 mm			183,100



Union



Union Elbow



Reducing Union



T-Type